



Recombinant Abnormal cell lineage protein 44 (lin-44)

Product Code	CSB-BP427112CXX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8XEH1
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	>85% (SDS-PAGE)
Sequence	SNDIPRTGG NKIVQPPKPN ILKQGCPSDL LHSRALRSIQ LACRSHPATV ISAFEGIQEG LQNCANLRF QQWDCSEAGN IMHDPPLLRQ GFRESSLIWA LSSASAAWGV ATACAQGWID DCACNNHMGQ NEYEFGGCTH GVQHGITASR KLLTKVGAVN SLLRKVEKHN LKAGRLAIKK TLISSCKCHG VSGSCQQKTC WKRTATLEHI TDYLVEKYAR AKLYTDDSVV KTTDLIYLEA SPDVCKVKS V AGRVCAWRNE THTQGDCDRL CCGNGFSIRH EVVRVKCDCE FVWCCNLVCK DCIQHRWIST CNGTPPKSLI F
Source	Baculovirus
Target Names	lin-44
Protein Names	Recommended name: Abnormal cell lineage protein 44 Alternative name(s): Wnt protein
Expression Region	22-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.