



# Recombinant Mouse 60S ribosomal protein L7 (Rpl7)

<b>Product Code</b>	CSB-BP020302MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P14148
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEAVPEK KKKK VATVPGTLKK KVPAGPKTLK KKVPAVPETL KKKRRNFAEL KVKRLRKKFA LKTLRKARRK LIYEKAKHYH KEYRQMYRTE IRMARMARKA GNFYVPAEPK LAFVIRIRGI NGVSPKVRKV LQLLRLRQIF NGTFVKLNKA SINMLRIVEP YIAWGYPNLK SVNELIYKRG YGKINKKRIA LTDNSLIARS LGKFGIICME DLIHEIYTVG KRFKEANNFL WPFKLSSPRG GMKKKTTTHFV EGGDAGNRED QINRLIRRMN
<b>Source</b>	Baculovirus
<b>Target Names</b>	Rpl7
<b>Protein Names</b>	Recommended name: 60S ribosomal protein L7
<b>Expression Region</b>	1-270
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L30P family of ribosomal proteins. It contains an N-terminal basic region-leucine zipper (BZIP)-like domain and the RNP consensus submotif RNP2. In vitro the BZIP-like domain mediates homodimerization and stable binding to DNA and RNA, with a preference for 28S rRNA and mRNA. The protein can inhibit cell-free translation of mRNAs, suggesting that it plays a regulatory role in the translation apparatus. It is located in the cytoplasm. The protein has been shown to be an autoantigen in patients with systemic autoimmune diseases, such as systemic lupus erythematosus. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.