



Recombinant Human Ubiquitin carboxyl-terminal hydrolase 17-like protein 5 (USP17L5)

Product Code	CSB-MP427555HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8MUK1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEDDSL YLRG EWQFNHFSKL TSSRPDAAFA EIQR TSLPEK SPLSCETRVD LCDDLAPVAR QLAPREKLPL SSRPAAVGA GLQNMGNTCY VNASLQCLTY TPPLANYMLS REHSQTCHRH KGCMLCTMQA HITRALHNPG HVIQPSQALA AGFHRGKQED AHEFLMFTVD AMKKA CLPGH KQVDHHSKDT TLIHQIFGGY WRSQIKCLHC HGISDTFDPY LDIALDIQAA QSVQQALEQL AKPEELNGEN AYHCGVCLQR APASKTLTLH TSAKVLILVL KRFSDVTGNK IAKNVQYPEC LDMQPYMSQP NTGPLVYVLY AVL VHAGWSC HNGHYFSYVK AQEGQWYKMD DAEVTASSIT SVLSQQAYVL FYIQKSEWER HSESVSRGRE PRALGAEDTD RRATQGELKR DHPCLQAPEL DEHLVERATQ ESTLDHWKFL QEQNKT KPEF NVRKVEG TLP PDVLVIHQSK YKCGMKNHHP EQQSSLLNLS SSTPTHQESM NTGTLASLRG RARRSKGKNK HSKRALLVCQ
Source	Mammalian cell
Target Names	USP17L5
Protein Names	Recommended name: Ubiquitin carboxyl-terminal hydrolase 17-like protein 5 EC= 3.4.19.12 Alternative name(s): Deubiquitinating enzyme 17-like protein 5 Ubiquitin thioesterase 17-like protein 5 Ubiquitin-specific-processing protease 1
Expression Region	1-530
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.