



Recombinant Arabidopsis thaliana SKP1-like protein 20 (ASK20)

Product Code	CSB-EP428490DOA-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8MQG7
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSEGLAVMK PETMKSYIWL QTADGSIQQV EQEVAMFCPM ICQEVIQKGV GSSKNHAISL PQRVNPAMFS LILDYCRFHQ LPGRSNKERK TYDERFIRMD TKRLCELTSA ADSLQLKPLV DLTSRALARI IEGKNPEEIR EIFHLPDDL EEEKLEPLKN SMDDPRI RLL NRYAKKRKE LKERERLKNV EVEEHVDERS VDDLFSFING RDHKAVKMSK GK KKKKKKKK QKIVSSNNIH DKESHDLRSK QQCVVEIGSS MREVPNLLSA EDDISTPNAG SEDEDIDDEI DPAMRELLDR EVEDFAQRLN SNWVRS LGKE RRPVHFSING NGTTRRHTGM VSIWRPCMFT LF
Source	E.coli
Target Names	ASK20
Protein Names	Recommended name: SKP1-like protein 20 Short name= AtSK20
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.