



Recombinant Neisseria meningitidis serogroup C Quinolinate synthase A (nadA)

Product Code	CSB-MP429373NEW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A9M2J9
Product Type	Recombinant Protein
Immunogen Species	Neisseria meningitidis serogroup C (strain 053442)
Purity	>85% (SDS-PAGE)
Sequence	<pre>MQTAARRSFD YDMPLIQTPT SACQIRQAWA KVADTPDRET ADRLKDEIKA LLKEKNAVLV AHYYVDPLIQ DLALETGGCV GDSLEMARFG AEHEAGTLVV AGVRFMGESA KILCPEKTVL MPDLEAEC SL DLGCP EEFAS AFCDQHPDRT VVVYANTSAA VKARADWVVT SSVALEIVSY LKSRGEKLIW GPDRHLGDYI RRETGADMLL WQGSCIVHNE FKGQELAALK AEHPDAVVLV HPESPQSVIE LGDVVGSTSK LLKAAVSRPE KKFIVATDLG ILHEMQKQAP DKEFIAAPTA GNGGCKSCA FCPWMAMNSL GGIKYALTSG RNEILLDRKL GEA AKLPLQR MLDFAAGLKR GDVFNGMGPA</pre>
Source	Mammalian cell
Target Names	nadA
Protein Names	Recommended name: Quinolinate synthase A EC= 2.5.1.72
Expression Region	1-370
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.