



Recombinant Mouse Annexin A6 (Anxa6)

Product Code	CSB-MP001847MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P14824
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	AKIAQGAMY RGSVHDFPEF DANQDAEALY TAMKGFSGSDK ESILELITSR SNKQRQEICQ NYKSPLYGKDL IEDLKYELTG KFERLIVNLM RPLAYCDAKE IKDAISGVGT DEKCLIEILA SRTNEQMHQL VAAYKDAYER DLESDIIGDT SGHFQKMLVV LLQGTRENDV VVSEDLVQQD VQDLYEAGEL KWGTDEAQFI YILGNRSKQH LRLVFDEYLK TTGKPIEASI RGELSGDFEK LMLAVVKCIR STPEYFAERL FKAMKGLGTR DNTLIRIMVS RSELDMLDIR EIFRTKYEKS LYSMIKNDTS GEYKALLKL CGGDDDAAGQ FFPEAAQVAY QMWELSAVSR VELKGTVCAA NDFNPDADAK ALRKAMKGIG TDEATIIDIV THRSNAQRQQ IRQTFKSHFG RDLMDLKSE ISGDLARLIL GLMMPPAHYD AKQLKKAMEG AGTDEKTLIE ILATRTNAEI RAINEAYKED YHKSLEDALS SDTSGHFRR LISLATGNRE EGGENRDQAQ EDAQVAEIL EIADTPSGDK TSLETRFMTV LCTRSYPHLR RVFQEFIKKT NYDIEHVIKK EMSGDVKDAF VAIVQSVKKN PLFFADKLYK SMKGAGTDEK TLTRVMVSR S EIDLLNIRRE FIEKYDKSLH QAIEGDTSGD FMKALLALCG GED
Source	Mammalian cell
Target Names	Anxa6
Protein Names	Recommended name: Annexin A6 Alternative name(s): 67 kDa calelectrin Annexin VI Annexin-6 Calphobindin-II Short name= CPB-II Chromobindin-20 Lipocortin VI Protein III p68 p70
Expression Region	2-673
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	Annexin VI belongs to a family of calcium-dependent membrane and phospholipid binding proteins. Although their functions are still not clearly defined, several members of the annexin family have been implicated in membrane-related events along exocytotic and endocytotic pathways. The annexin VI gene is approximately 60 kbp long and contains 26 exons. It encodes a protein of about 68 kDa that consists of eight 68-amino acid repeats separated by linking sequences of variable lengths. It is highly similar to human annexins I and II sequences, each of which contain four such repeats. Exon 21 of annexin VI is alternatively spliced, giving rise to two isoforms that differ by a 6-amino acid insertion at the start of the seventh repeat. Annexin VI has been



implicated in mediating the endosome aggregation and vesicle fusion in secreting epithelia during exocytosis.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.