



# Recombinant Methanococcus maripaludis Signal recognition 54 kDa protein (srp54)

<b>Product Code</b>	CSB-BP432530MNQ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	A9A9B0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Methanococcus maripaludis (strain C6 / ATCC BAA-1332)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLDKLGQNLSEALNKIKNATFVDKKLVEV IKDIQKALIQSDVNVKLVFN MSKEIERKAI EEAPPKGLSK KEHIVKIVYD ELVKLLGETI QKLELDASKK SVILLIGIQG SGKTTSAAKL ARYIQKKGLR PGLIAADVYR PAAYQQLKQL SEKINVPLFG DETRTKTPVD ITKEGMEKLL KVDVIIDTA GRHKEEEGLL AEMKEMKDLT NPNEILVID GTLGQQAKNQ AKAFKDSVSE IGSILVTKLD GSAKGGGALS AVAEINAPIK FIGTGEGVDN LEQFDPKKFI SRILGFGDLD SLLEKTEDIM DESAEESIDS ILKKGFTLIE LYAQLETISK MGPMKQILSM IPGMGGNMPK EAAQLTEDKL KRYKIMMDSM TMEKENPEL IKTSRLQRIA RGAGVKPEEI KDLLKYSTT KNAFGNLKRG KMPKMGGQMG QIMRQLMYKD
<b>Source</b>	Baculovirus
<b>Target Names</b>	srp54
<b>Protein Names</b>	Recommended name: Signal recognition 54 kDa protein Alternative name(s): SRP54
<b>Expression Region</b>	1-450
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.