



Recombinant Methanococcus maripaludis Signal recognition 54 kDa protein (srp54)

Product Code	CSB-YP432530MNQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A9A9B0
Product Type	Recombinant Protein
Immunogen Species	Methanococcus maripaludis (strain C6 / ATCC BAA-1332)
Purity	≥85% (SDS-PAGE)
Sequence	MLDKLGQNLSEALNKIKNATFVDKKLVEV IKDIQKALIQSDVNVKLVFN MSKEIERKAI EEAPPKGLSK KEHIVKIVYD ELVKLLGETI QKLELDASKK SVILLIGIQG SGKTTSAAKL ARYIQKKGLR PGLIAADVYR PAAYQQLKQL SEKINVPLFG DETRTKTPVD ITKEGMEKLL KVDVIIDTA GRHKEEEGLL AEMKEMKDLT NPNEILVID GTLGQQAKNQ AKAFKDSVSE IGSILVTKLD GSAKGGGALS AVAEINAPIK FIGTGEGVDN LEQFDPKKFI SRILGFGDLD SLLEKTEDIM DESAEESIDS ILKKGFTLIE LYAQLETISK MGPMKQILSM IPGMGGNMPK EAAQLTEDKL KRYKIMMDSM TMEEKENPEL IKTSRLQRIA RGAGVKPEEI KDLLKYSTT KNAFGNLKRG KMPKMGGQMG QIMRQLMYKD
Source	Yeast
Target Names	srp54
Protein Names	Recommended name: Signal recognition 54 kDa protein Alternative name(s): SRP54
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.