



Recombinant Human Pregnancy-specific beta-1-glycoprotein 1 (PSG1)

Product Code	CSB-EP018849HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P11464
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>QVTIEA EPTKVSEGKD VLLL VHNL PQ NLTGYI WYKG QMRDLYHYIT SYVVDGEIII YGPAYSGRET AYSNASLLIQ NVTREDAGSY TLHIIKGDDG TRGVTGRFTF TLHLETPKPS ISSSNLNPRE TMEAVSLTCD PETPDASYLW WMNGQSLPMT HSLKLSETNR TLFLLGVTKY TAGPYECEIR NPVSASRSDP VTLNLLPKLP KPYITINLN PRENKDVLNF TCEPKSENYT YIWWLNGQSL PVSPRVKRPI ENRILILPSV TRNETGPYQC EIRDYGGIR SDPVTNLNLY GPDLPRIYPS FTYYRSGEVL YLSCSADSNP PAQYSWTINE KFQLPGQKLF IRHITTKHSG LYVCSVRNSA TGKESSKSMT VEVSDWTV P</p>
Source	E.coli
Target Names	PSG1
Protein Names	Recommended name: Pregnancy-specific beta-1-glycoprotein 1 Short name= PS-beta-G-1 Short name= PSBG-1 Short name= Pregnancy-specific glycoprotein 1 Alternative name(s): CD66 antigen-like family member F Fetal liver non-spe
Expression Region	35-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.