



Recombinant Human Complement C1q and tumor necrosis factor-related protein 9B (C1QTNF9B)

Product Code	CSB-MP450719HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B2RNN3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	Q DTCRQGHPIG PGNPGHNGLP GRDGRDGAKG DKGDAGEPGC PGSPGKDGTS GEKGERGADG KVEAKGIKGD QGSRGSPGKH GPKGLAGPMG EKGLRGETGP QQQKGNKGDV GPTGPEGPRG NIGPLGPTGL PGPMGPIGKP GPKGEAGPTG PQGEPGVRGI RGWKGDRGEK GKIGETLVLP KSAFTVGLTV LSKFPSSDVP IKFDKILYNE FNHYDTAVGK FTCHIAGVYY FTYHITVFSR NVQVSLVKNG VKILHTRDAY VSSEDQASGS IVLQLKLGDE MWLQVTGGER FNGLFADEDD DTTFTGFLLF SSQ
Source	Mammalian cell
Target Names	C1QTNF9B
Protein Names	Recommended name: Complement C1q and tumor necrosis factor-related protein 9B Alternative name(s): C1q/TNF-related protein 9B Short name=CTRP9B Complement C1q and tumor necrosis factor-related protein 9-like
Expression Region	20-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.