



Recombinant *Saccharomyces cerevisiae* tRNA (guanine-N (7)-)-methyltransferase subunit TRM82

Product Code	CSB-YP451433SVP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3LGB9
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSVIHPLQNL LTSRDGSLVF AIIKNCILSF KYQSPNHWEF AGKWSDDDFDK IQESRNTTAK EQQGQSSENE NENKKLKS NK GDSIKRTAAK VPSPGLGAPP IYSYIRNLRL TSDESRLIAC ADSDKSLLVF DVDKTSKNVL KLRKRFCFSK RPNASIAED DTTVIIADKF GDVYSIDINS IPEEKFTQEP ILGHVSMLTD VHLIKDSDGH QFIITSDRDE HIKISHYPQC FIVDKWFLFGH KHFVSSICCG KDYLLLSAGG DDKIFAWDWK TGKNLSTFDY SSLIKPYLND QHLAPPRFQN ENNDIIEFAV SKIIKSKNLP FVAFFVEATK CIIILEMSEK QKGDALALKQI ITFPYNVISL SAHNDEFQVT LDNKESGQVQ KNFAKFIEYN LNENSFVVNN EKSNEFDSAI IQSVQGDSNL VTKKEEIYPL YNVSSLRKHG EHYS
Source	Yeast
Target Names	TRM82
Protein Names	Recommended name: tRNA (guanine-N(7)-)-methyltransferase subunit TRM82 Alternative name(s): Transfer RNA methyltransferase 82
Expression Region	1-444
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.