



Recombinant *Saccharomyces cerevisiae* Spindle pole component 29 (SPC29)

Product Code	CSB-EP451440SVP-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3LKV0
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MDYSNFGNSA SKKFQDDTLN RVRKEHEEAL KKLREENFSS NTSELGNKKH YRAQERMSSP LHRLSPTGKS DDRKVKSPLD DKLRRQLREG NTRLPPPPFS SYGMPPTNRS NLDRI RRRTS SPVRTDKFAS QNVIDDQRLE IKYLERIVYD QGTVIDNLTS RITRLESFIL NSISDRGDKN FASLEHSRSF SGFPTNKTYG LQMGGLYEND MPYRRSSDNI NKEGAREDRS SQIHIENEST EDILKILSSS FHN
Source	E.coli
Target Names	SPC29
Protein Names	Recommended name: Spindle pole component 29
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.