



Recombinant *Drosophila ananassae* Queuine tRNA-ribosyltransferase subunit QTRTD1 homolog (GF24662)

Product Code	CSB-MP451460DKY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3MA91
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MKLVIESISK NSGRLGHLRI REGGPDLKTP LLLQTTKGGG IPWLAGDVFE HQVSEQEPQVL ELTSLNMGHM ADALQQWNNQ GRGVSDYVGF PGNPNLLLLR DPCETTPSGS NDRDIQPLFT RRGRESLTAE TYIEMVAKIQ PDIFQGLSDA DTNAESSKKR AQKSVDRTK FMQYIYEQSG KVRGTLAPI VGGYNTFART QSIKHAREQP AESYGGYVLE GFHNGLSAT SLDTSKLLPI VEHCVKQLEE DKPRIIPGAY TPLTVLELIR LGIDVFDTSY AYCATHNFKA LTFTFVKNEK EDKPFLDITD DTIKEEFTPL LRDCSCLTCQ KHTRAYLHHL YKTNELLGPI LLMIHNLHHY MAFFERIREA VANDTLPALT ELVKTQNGKS SSDFSIVQNT KVISKATMSK GFVAAAV
Source	Mammalian cell
Target Names	GF24662
Protein Names	Recommended name: Queuine tRNA-ribosyltransferase subunit QTRTD1 homolog EC= 2.4.2.29 Alternative name(s): Queuine tRNA-ribosyltransferase domain-containing protein 1
Expression Region	1-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.