



Recombinant Zea mays Ninja-family protein 7

Product Code	CSB-BP451826ZAX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4FXN6
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MDDDNGLLELS LGLSLGGSSG KAKARDAPLE PKAEPQVEES SSKGGLQTPD APSGKFYQKS ADNHEQNSKQ RHSPVAPQFG NFWGQPGNSS GPVVDGSVEP VSHQPQLSSY QDGRMPINSG NNSEEQKPVS GNHKL PSEEM NFQKKHQTAS DQPDAFSKSS DGGAKNAPIS ISTDDGSTGE NEDVAESEAE GSNSWLVAQR EDSAKGSVVN KGTDKRKRSAD AAADGFKGKR QPSFSGSESS SGKLTGPNLI SMQASNVAL PYQVPGQVSG PPVLTNAPNF HPVCTVQLRP PTNGGLAVQT MGSASQVAFG YPAVQLPTLE TGSSWAFGAP PQALSSFTAK DKADQTGTKQ VDDGKKPREA GASSSAHAED EKKADRGLSL MGSAIRPGIA PNVKFGGSGS YPDLPWVSTI GAGPNGRTIS GVTYKFRNE VKIVCACHGT HMSPEEFMRH ASADAPAQDN SATLPAFPAG NQATSAEN
Source	Baculovirus
Protein Names	Recommended name: Ninja-family protein 7
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.