



Recombinant *Drosophila yakuba* Cytosolic Fe-S cluster assembly factor NUBP1 homolog (GE12759)

Product Code	CSB-EP451897DMR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4P9A8
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila yakuba</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MQAPPPEHCP GVESEQAGKG SACAGCPNQG VCSDPNKKLE DPGKALVVES MKDVSNKLLI LSGKGGVVGKS TVTSLLLTRYL ARSNPDSNFG VLDIDICGPS QPRLMGALGE SVHQSGYGWS PVGIEDNVCL MSIGFLLDSV DDAIWRGPK KNGMIRQFLS EVDWGNLDDL LLDTPPGTSD EHLSVVSYLK DDTNPESLRA VMVTTPQEVS LLDVRKEINF CKKQNIPIVG VIENMSSFRC GHCNNSSEIF PAKTGGAAAM CAEMGIPLL G SLPLDQQISK ACDSGEDLTE FKNVTTEALE GICSKIMASF S
Source	E.coli
Target Names	GE12759
Protein Names	Recommended name: Cytosolic Fe-S cluster assembly factor NUBP1 homolog
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.