



Recombinant *Drosophila yakuba* Molybdopterin synthase sulfur carrier subunit (Mocs2)

Product Code	CSB-EP451903DMR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4PUD0
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila yakuba</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MNADGPVVNV HVLFFAKSRE LANTPRSKVD LPTEITASDL LELLVSKFGL TPIRDNLILA HNESYIDNLS ERILFKEGDE LAIIPPLSGG
Source	E.coli
Target Names	Mocs2
Protein Names	Recommended name: Molybdopterin synthase sulfur carrier subunit Alternative name(s): Molybdenum cofactor synthesis protein 2 small subunit Molybdenum cofactor synthesis protein 2A Short name= MOCS2A Sulfur carrier protein MOCS2A
Expression Region	1-90
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.