



# Recombinant Rat Calcium/calmodulin-dependent protein kinase type II subunit gamma (Camk2g)

<b>Product Code</b>	CSB-BP004468RA
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P11730
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MATTATCTRF TDDYQLFEEL GKGAFSVVRR CVKKTSTQEY AAKIINTKKL SARDHQKLER EARICRLKLN PNIVRLHDSI SEEGFHYLVF DLVTGGELFE DIVAREYYSE ADASHCIHQI LESVNHIIHQH DIVHRDLKPE NLLLASKCKG AAVKLADFGL AIEVQGEQQA WFGFAGTPGY LSPEVLRKDP YGKPVDIWAC GVILYILLVG YPPFWDEDQH KLYQQIKAGA YDFPSPEWDT VTPEAKNLIN QMLTINPAKR ITADQALKHP WVCQRSTVAS MMHRQETVEC LRKFNARRKL KGAILTTMLV SRNFSVGRQS SAPASPAASA AGLAGQAACS LLNKKSDGGV KKRKSSSSSVH LMPEQTTVVH NATDGIKGST ESCNTTTEDE DLKVRKQEII KITEQLIEAI NNGDFEAYTK ICDPGLTSFE PEALGNLVEG MDFHKFYFEN LLSKNSKPIH TTILNPHVHV IGEDAACIAY IRLTQYIDGQ GRPRTSQSEE TRVWHRRDGGK WLVNHYHCSG APAAPLQ
<b>Source</b>	Baculovirus
<b>Target Names</b>	Camk2g
<b>Protein Names</b>	Recommended name: Calcium/calmodulin-dependent protein kinase type II subunit gamma Short name= CaM kinase II subunit gamma Short name= CaMK-II subunit gamma EC= 2.7.11.17
<b>Expression Region</b>	1-527
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	The product of this gene belongs to the Serine/Threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. Calcium signaling is crucial for several aspects of plasticity at glutamatergic synapses. In mammalian cells the enzyme is composed of four different chains: alpha, beta, gamma, and delta. The product of this gene is a gamma chain. Six alternatively spliced variants that encode six different isoforms have been characterized to date. Additional alternative splice variants that encode different isoforms have been described, but their full-length nature has not been determined.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.