



# Recombinant Pig Neurofilament heavy polypeptide (NEFH)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP015686PI-B                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Uniprot No.</b>       | P12037                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Sus scrofa (Pig)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | MRGAVLRLGA ARGQLRLEQE HLLEDIAHVR QRLDDEARQR QEAEAAARAL<br>ARFAQEAEAA RVELQKKAQA LQEECGYLRR HHQEEAQAEA RDALKCDVTS<br>ALREIRAQLE GHAVQSTLQQ EEWFRVRLDR LSEAAKVNTD AM                                                                                                                                                                                                                                                                                                                                          |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Target Names</b>      | NEFH                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Protein Names</b>     | Recommended name: Neurofilament heavy polypeptide Short name= NF-H<br>Alternative name(s): 200 kDa neurofilament protein Neurofilament triplet H protein                                                                                                                                                                                                                                                                                                                                                    |
| <b>Expression Region</b> | 1-142                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | Full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Target Details</b>    | Neurofilaments are type IV intermediate filament heteropolymers composed of light, medium, and heavy chains. Neurofilaments comprise the axoskeleton and functionally maintain neuronal caliber. They may also play a role in intracellular transport to axons and dendrites. This gene encodes the heavy neurofilament protein. This protein is commonly used as a biomarker of neuronal damage and susceptibility to amyotrophic lateral sclerosis (ALS) has been associated with mutations in this gene. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                 |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                          |