



Recombinant Danio rerio Protein TILB homolog (Irrc6)

Product Code	CSB-EP453211DIL-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3DH20
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MVRISEDLR RRAEHNNGEI FSLEELSLHQ QDIQRIEHIH KWCRDLKILY LQNNLIPKIE NVGRLKKLEY LNLALNNIEV IENLEGCESL QKLDLTVNFV GRLSSVETLK HNLHLKELYL VGNPCAHEYQG YRQYVVATVP QLQSLDGKEI SRAERIQALQ ELDAVRTRVL QQETKYLEER EKQKSNANEH PEINQSLSES QNGTQQYPES SSKTHTEAED EEREFWEKPC PFTPESRLEA HRHLEEKRRA NEKEKEKPKT KTPRTLITPD GRVLNVNEPK LDFSLFEDEN NCLLLDLHVV RHMDSSLLDV DVQPMYVRVT VKGKVFQLVL PAEVKPDSSS AQRSQTGHL LLILPLANED VKPKKRTIRP TSVTSNQNNK KDTRAAPRRE LLEVDPGLAG SLANIVPKGQ ESSHNPQRCG LEERPVSDF VDDPEVPLM
Source	E.coli
Target Names	Irrc6
Protein Names	Recommended name: Protein TILB homolog Alternative name(s): Leucine-rich repeat-containing 6-like protein Leucine-rich repeat-containing protein 6 Seahorse
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.