



Recombinant Danio rerio Protein TILB homolog (Irrc6)

| | |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-YP453211DIL |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B3DH20 |
| Product Type | Recombinant Protein |
| Immunogen Species | Danio rerio (Zebrafish) (Brachydanio rerio) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MVRISEDLR RRAEHNNGEI FSLEELSLHQ QDIQRIEHIH KWCRDLKILY LQNNLIPKIE NVGRLKKLEY LNLALNNIEV IENLEGCESL QKLDLTVNFV GRLSSVETLK HNLHLKELYL VGNPCAHEYQG YRQYVVATVP QLQSLDGKEI SRAERIQALQ ELDAVRTRVL QQETKYLEER EKQKSNANEH PEINQSLSES QNGTQQYPES SSKTHTEAED EEREFWEKPC PFTPESRLEA HRHLEEKRRA NEKEKEKPKT KTPRTLITPD GRVLNVNEPK LDFSLFEDEN NCLLLDLHVV RHMDSSLLDV DVQPMYVRVT VKGKVFQLVL PAEVKPDSSS AQRSQTTGHL LLILPLANED VKPKKRTIRP TSVTSNQNNK KDTRAAPRRE LLEVDPGLAG SLANIVPKGQ ESSHNPQRCG LEERPVSDF VDDPEVPLM |
| Source | Yeast |
| Target Names | Irrc6 |
| Protein Names | Recommended name: Protein TILB homolog Alternative name(s): Leucine-rich repeat-containing 6-like protein Leucine-rich repeat-containing protein 6 Seahorse |
| Expression Region | 1-440 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |