



Recombinant *Saccharomyces cerevisiae* ASTRA-associated protein 1 (ASA1)

Product Code	CSB-MP453388SVP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3LK55
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MRGFSNEIIL KRTLTLSDFT LRYHKGITA LQVIKAPSVS NVPVLLSGDN YGYFVMWDLV TKRPITHIEI EGNSHIAFW WVETTNVLYI LSKDSMLRIF ELDSSTQLSI DLVRKLSQAN KTDHLQWTKI YEMPINTLNF ANFIIAEVK PTKDNKSYRL VCCHTDDSET IDIYQIIEDS TFKLKRPFNN INFPRFLKQQ NFLGISKDSK FGIIMRFAKL NDVIFLGYEN GFVVGFKITF DEGLQRDIAE LVHVSNDHYP NPILDMCVSG DELYSCSTDD FITKYKIPVN LQLETKYLRD DALLIKCPSS LRVSEPSKVH LPLKNIGHID KVKDDYLVVS SWSGMTIVYN MRTSEVEQTF VKSKNNLVVS DSSMGDLTNG SGSNTESSSK SHNYKVGAMT CLESFDVQSD GLRLGQLRRI KALAKCNWCL IGYEDGTIKL NKI
Source	Mammalian cell
Target Names	ASA1
Protein Names	Recommended name: ASTRA-associated protein 1
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.