



# Recombinant *Saccharomyces cerevisiae* Protein SWT21 (SWT21)

<b>Product Code</b>	CSB-BP453394SVP
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B3LP36
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEKKVICQDI FWSCDGTSFV SVHNDFGIRQ YLVPEESNTD KLNRLNLLLPF TRFFRNQSIV SCAIDPFYTL YNENSDRLAG DRIVVGGKNF PLQLYSLMDG QCILSYDTMN KINGEYETVY SVKIDVESRV YTGSCRNKVA IYDKSRRDAV WMNQSTKKAS KGRQSIISCF EEQPMGGQAL SRGSLLCGSY ANEMFQVDCR HQRLERLNYT RTVAGGIVQI LTSDNGRYVY VVRRNSDAIS IYDRRLQHE LNLVRLPFRI HHNSAKLKAY IDTAYGLSMG TPQGILNWG RDLVEFGGVP SHNSVEDPLI TSIPPESEWR TNLDSTIPAT VVKNCPGDPE LFALSHGGTI SLCRFGG
<b>Source</b>	Baculovirus
<b>Target Names</b>	SWT21
<b>Protein Names</b>	Recommended name: Protein SWT21 Alternative name(s): Synthetic With TGS1 protein 21
<b>Expression Region</b>	1-357
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.