



Recombinant *Drosophila ananassae* Structure-specific endonuclease subunit SLX1 homolog (slx1)

Product Code	CSB-BP453408DKY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3M0F3
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSFRKFTATA DPEKDETIAR KGHFYGVYLL CSQSLDSRYR GKCYVGFTVN PKRRIKQHNH R GCDFFGAKKT SRKGPWQMVM IVHGFPNNIV ALQFEWAWQQ PTLSTRLKIF PELKRKNPRE SHFDYNFRIL NRMLGVGPWN RLALKIRWLE TDYERGFVLP LPRHMEIVSG KVSISSSQRR KGEDAAATVP VAWAPECHLC MQRIDQPERS RLGCTNPTCR LTCHMLCLAN YLLGDEPGHY IPVGGECPLC ETRLSWSALL QRKRLNGIP EELLDQEEAD EDLSDGPDVD SDVEVLESD
Source	Baculovirus
Target Names	slx1
Protein Names	Recommended name: Structure-specific endonuclease subunit SLX1 homolog EC= 3.6.1.-
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.