



Recombinant Danio rerio Polycomb group RING finger protein 5-A (pcgf5a)

Product Code	CSB-EP455182DIL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3DK16
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MATHRKHLVR DFNRYITCSI CRGYLIKPTA VTECLHTFCK SCIVQHFEES NECPECGIQV HETNPLEMLR LDKTLEEIIF KLVPLREKE EHQESEFWRK HKIKSNGEDG PRAKKSRLSG EDDDGNGGDY HRSDPQIAIC LDCLRNNGQS GESIVKGLMK KFIRCSTRVT VGTIKKFLCV KLKLPSSYEL DVLCNGEIMG KDHTLEFIYR TRWRLQGDSA YPMVLEYRPR IDFG
Source	E.coli
Target Names	pcgf5a
Protein Names	Recommended name: Polycomb group RING finger protein 5-A
Expression Region	1-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.