



Recombinant *Drosophila ananassae* WD repeat-containing protein 55 homolog (GF18399)

Product Code	CSB-BP455379DKY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3M1G0
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MHTHNHFKTP SDAEEVDDLD DEMVMGVIAE IEQEVLFESD SDDDGFGAEQ LGAPPPDDDN DQISSSDDSF DPNAEDSDSD DSMQHEEQAQ ATSAKRRKDD DGPSGSNRDV ETNFDLDLED ETDETVRAII AAIKKPRSAP PEIKLEDFVT DISFHPERNI IALATIIGDV HLYEYANEGN KLIRTIEVHS KACRDVEFTE DGRNLLTCSK DKCVMVTDME TEKLKLYET AHDDAINTLH VLNENLFATG DDAGTVKLWD LRTKQHVFE LKQIDDQVTQL LSNEQNTLLL ATSADGYLTT FNIPGRKLYV QSEPYEEELN CMGIYRGDSK LVVGTSGKGL YSYNWGSFGY HCDMYPGIKS PISLMIPITD RIACVAGEDG NIRACHITPY RNLGVVGQHN MPIEALDVNS TGELLASSSH NNDVRFWNVK YFEDFGDIKY NEKHNAYKEQ RHNLPSSKCT NTGDIFFADLA KQEDDEDDDA T
Source	Baculovirus
Target Names	GF18399
Protein Names	Recommended name: WD repeat-containing protein 55 homolog
Expression Region	1-481
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.