



Recombinant *Drosophila erecta* Cytosolic Fe-S cluster assembly factor NUBP1 homolog (GG22765)

Product Code	CSB-MP455396DLK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3NNJ9
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MQAPPPEHCP GLESEQAGKA SACAGCPNQG VCSDPNKKLE DPGKALVVES MKDVKHKLLI LSGKGGVGKS TVTSLTRYL ARSNPDSNFG VLDIDICGPS QPRLMGALGE SVHQSGYGWS PVGIEDNVCL MSIGFLLGSV DDAIWRGPK KNGMIRQFLS EVDWGNLDDL LLDTPPGTSD EHLSVVSYLK DDTNPELCA VMVTTPQEVS LLDVRKEINF CKKQNIPIVG VIENMSSFRC GNCGNSSSEIF PAKTGGAAAM CAEMGIPLLG SLPLDQQIAK ACDSGEDITE FKNVTTEALD GICSKIIASF S
Source	Mammalian cell
Target Names	GG22765
Protein Names	Recommended name: Cytosolic Fe-S cluster assembly factor NUBP1 homolog
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.