



Recombinant Zea mays Ninja-family protein 1

Product Code	CSB-EP455717ZAX-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4FAF3
Product Type	Recombinant Protein
Immunogen Species	Zea mays (Maize)
Purity	>85% (SDS-PAGE)
Sequence	MEGFSRDLLC GIGKGDAPPP EKRPGLREE MEEVELSLGL SLGGRFGLDR KGSKLPRSSS VAAMLTTPVE VPAPPSLSRA SSLPVQAEAS EVERKQGLDG WGSCREGGGL GVQHAARLPA SGNPSSASSV GEGQILQGTL MRTSSLPAVI EASGNDDWKK RKEAQLKRL EVKKKRIERR NSLACNTSKE AAGQSPKEMN ANTDKLVSSD ETIVSANESH SSGKHLVKGL PPKYQATITS EDSSSAMRKK PNSAFKGTAI TEEQNSSSSV PSSGEAISSV TAPSLPLLSL VPITATLGSR EDQSILGRAG ARANGMGDVE RRMMQEMPGV FTKGLSNGSR VEGLYKYSK GEVRIVCICH GSFLTPSEFV EHAGAGKVDN PLRHIVVSAT PNL
Source	E.coli
Protein Names	Recommended name: Ninja-family protein 1
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.