



Recombinant *Drosophila sechellia* Ribosome biogenesis protein WDR12 homolog (GM11446)

Product Code	CSB-EP455734DMG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4HWV6
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila sechellia</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDVDNGEGQV QVHLKTKQEQ YAVPDVPYAI DGTVTTVELN TFDVALLRQK DGSSDTEFDF LVFDEYLRGR LCDHLREKAI SFEDAIEIEY VERFPAPEPQ DCLLHDDWVS AVKARGKWIL SGCYDNSLNL WTNKKGKHILT ISGHTAPIKA VDWISLDEET GRFVSTSQDQ TAMLWKWNVG SNAVDCVSVC KHERGVDSV SVSPDGLRFA TGSWDTMLKV WSAELDDGVE GSSKRMKESG VRTPKITLQG HRESVSAVQW MDATLLTGS WDYTLKVWDL SLEGIKTEIS TNKSIFDASY SKLNRLILTA SADKNLRLYD PRTNQGSVVR NTYLGHNWV QTVMWSTTEE FLFVSGAYDN QNKLWDCRSP KAPLYDLLGH GDKVLDIDWS NPKYIVSGGV DNTVRVFKSR KALAEDETETK
Source	E.coli
Target Names	GM11446
Protein Names	Recommended name: Ribosome biogenesis protein WDR12 homolog
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.