



Recombinant *Drosophila sechellia* Structure-specific endonuclease subunit SLX1 homolog (slx1)

Product Code	CSB-MP455736DMG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4I3R2
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila sechellia</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MNSYDPQDTG SQQEETVALK GHFYGVYLLC SQSLDPRYRG KCYVGFTVNP KRRIRQHNLG CDFGGARKTS RKG PWLMVMI VHGFNNNTVA LQFEWAWQQP SLSTR LKMYP ELKRKLPRET FFDYNFRILS HMLGVGPWNR LPLTVRWLET DYERPFSKTL PNHMEIVSGK VSISASQRRR PDDAVPPPPV AWASECHLCM QQIEQPEKSR LGCTNQMCRL TCHMVCLANY LLGDETGHYI PVGGECP LCE TRLSWAALLQ RKRLLMGVPE ELQDHD E DLS DDIDV DSDIE DGAQELSN
Source	Mammalian cell
Target Names	slx1
Protein Names	Recommended name: Structure-specific endonuclease subunit SLX1 homolog EC= 3.6.1.-
Expression Region	1-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.