



Recombinant Drosophila sechellia Protein king tubby (king-tubby)

Product Code	CSB-BP455737DMG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4I7J9
Product Type	Recombinant Protein
Immunogen Species	Drosophila sechellia (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MEAYIRQKRA SPGMVQASDL QINRPMSGMR SNSRELHAYD GPMQFISSPQ NPDQILTNGS PGGINSVAMN TSRNHSNNMR SLSTINQEAD LIEEISSHEL EDEESSPVTV IEQQQQSASH SANSTQSQKP RARQHSFSDN LDEDDYTNRN VAGAAPVRPA GMASSPYKDA TLDGSSNGTG NGTGGESEGD VIGNIDQFVM QPAPQGVLYK CRITRDRKGM DRGLFPIYYL HLERDYGKKI FLLGGRKRKK SKTSNYIVSC DPTDLARNAD GFCGKLRNV FGTSFTVFDN GNKESTESPR LDLAVIYDT NILGFKGPRN MTVILPGMTE DDQRVKISSA DPKQQGILD WKMKNMNDNIV ELHNKTPVWN DETQSYVLNF HGRVTQASVK NFQLVHSDP EYIVMQFGRT SEDVFTMDYR YPLCAMQAFA IALSSFDGKI ACE
Source	Baculovirus
Target Names	king-tubby
Protein Names	Recommended name: Protein king tubby
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.