



Recombinant Drosophila yakuba Ribosome biogenesis protein WDR12 homolog (GE13328)

Product Code	CSB-BP455782DMR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4P116
Product Type	Recombinant Protein
Immunogen Species	Drosophila yakuba (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDVENGEGQV QVHLKTKQEH YAVPDVPYAI DGTVTTVELN TFDVALLRQK DGSSDTEFDF LVFDEYLRGR LCDHLREKAI SFEDAIEIEY VERFPAPEPQ DCLLHDDWVS AVKASGKWIL SGCYDNTLNL WTNKGKHILT ISGHTAPIKA VDWISLDEET GRFVSTSQDQ TAMLWQWNVG SNSVECVSVC KGERGVDSV SVSPDGLRFA TGSWDTMLKV WSAEVEDAVE GSSKRMKESG VRTPKITLQG HRESVSAVQW MDASTLLTGS WDHTLKVWDL SLEGIKTEIS TNKSIFDASY SKLNRLILTA SADKNLRLYD PRTNQGSVVR NTYLGHNWV QTVMWSTTEE FLFVSGAYDN QNKLWDCRSP KAPLYDLLGH GEKVLIDIDWS NPKYIVSGGV DNTVRVFKSR KALADDAETK
Source	Baculovirus
Target Names	GE13328
Protein Names	Recommended name: Ribosome biogenesis protein WDR12 homolog
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.