



Recombinant *Drosophila simulans* tRNA (guanine-N (7)-)-methyltransferase

Product Code	CSB-EP455805DMJ-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4R338
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila simulans</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MVATGGQAQD QSQNQEPDVL NPTSAVTGLP QKRYRQRAH SNPIADHSFD YPARPEDVDW RSMYPSIQQG QQVSFADVGC GYGGFLVTLG EMFPEKLSIG MEIRVKVSDY VVDRIAALRR RSADTGAYQN IACLRTNAMK YLPNYFSKGQ LEKMFFLYPD PHFKRAKHKW RIINQALLSE YAYVLRKGGGL VYTM TDVEDL HKWIVTHMEE HPLYERLTEE EANADPITPK LYQSSEEGAK VVRNKGDHFL AIFRRL
Source	<i>E.coli</i>
Target Names	GD16394
Protein Names	Recommended name: tRNA (guanine-N(7)-)-methyltransferase EC= 2.1.1.33 Alternative name(s): tRNA (guanine(46)-N(7))-methyltransferase tRNA(m7G46)-methyltransferase
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.