



Recombinant *Saccharomyces cerevisiae* Mitochondrial zinc maintenance protein 1, mitochondrial (MZM1)

Product Code	CSB-EP457256SVP-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3LFH4
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	DAEVLL AARAKMRSGM LCPPDPKLT EDQIQHLEDV AVFLRRNLVQ GKKVDGSSTK EPRYHLNIHK DTELGDNETI ADPTARVKTN LKARPFKCS D KKQ
Source	E.coli
Target Names	MZM1
Protein Names	Recommended name: Mitochondrial zinc maintenance protein 1, mitochondrial Alternative name(s): Altered inheritance of mitochondria protein 8 Found in mitochondrial proteome protein 36
Expression Region	25-123
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.