



Recombinant *Drosophila ananassae* Ubiquitin carboxyl-terminal hydrolase calypso (calypso)

Product Code	CSB-BP457289DKY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3MIV9
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MNMTGGGAGP QAGPNSSCNI PNSSLATAP VAATMPIAQL ADGWLELESD PGLFTLLLED FGCHDVQVEE VYDLQKPIES PYGFIFLFRW IEERRARRKI VETTAEIFVK DEEAISSIFF AQQVVPNSCA THALLSVLLN CNENNLQLGD TL SRLKAHTK GMSPENKGLA IGNTPELACA HNSHAMPQAR RRLERTGAGV SSCRFTGEAF HFVSFVPING QLFELDGLKP YPMNHGGWED HEDWTDKFRR VMTERLGIAT GEQDIRFNLM AVVPDRRIAI THKLKMLRTN QAI VSGTLQK LLKADEQGEG GNGDPQR PDT PSTLLEPSAF TARDLQSLK NLDTEIAINE QHLADENDRR HKFKVDASRR THNYDKFICT FLSMLAHQGV LGELVSQHLL PSKKISGQSA ANRLNKQ SNA TANAGATAAG AAGAAPKSQQ QQAAAAKNGK SPSKTPGRRR KGRNKCKKRK
Source	Baculovirus
Target Names	calypso
Protein Names	Recommended name: Ubiquitin carboxyl-terminal hydrolase calypso EC=3.4.19.12 Alternative name(s): BAP1 homolog
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.