



# Recombinant Drosophila simulans Eukaryotic translation initiation factor 3 subunit G-2 (eIF3-S4-2)

<b>Product Code</b>	CSB-EP457722DMJ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B4QS37
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Drosophila simulans (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKSFITSWAD EVDADYVDGL PPSNEYIKGD YKYVTEYKFN DDGKKVKVVR TFKIEKQIVP KAVARRRNWV KFGDSRSDKP GPNSQTTMAS EEILMQFIGS KEFDQTHETQ LDPGKNIACR RICNGEHWSV NCPYKGTSMR SKTMMETKAN AAAAAISDP SKTGKYVPPF MKDGGGSGS KNWGRGRDRD DSSAVRISNL SESMTETDLE ELVKKIGPHT KMYLAREKNS GLCKGFAYVH FKFRQDAAAA IEVLNGHGYD HLILCVEWSK PQP
<b>Source</b>	E.coli
<b>Target Names</b>	eIF3gb
<b>Protein Names</b>	Recommended name: Eukaryotic translation initiation factor 3 subunit G-2 Short name= eIF3g-2 Alternative name(s): Eukaryotic translation initiation factor 3 RNA-binding subunit 2 Short name= eIF-3 RNA-binding subunit 2 Eukaryotic
<b>Expression Region</b>	1-273
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.