



Recombinant Human Zinc finger protein 35 (ZNF35)

Product Code	CSB-MP026693HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P13682
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTAELREAMA LAPWGPVKVK KEEEEENFP GQASSQQVHS ENIKVWAPVQ GLQTGLDGSE EEEKGQNISW DMAVVLKATQ EAPAASTLGS YSLPGTLAKS EILETHGTMN FLGAETKNLQ LLVPKTEICE EAEKPLIISE RIQKADPQGP ELGEACEKGN MLKRQRIKRE KKDFRQVIVN DCHLPESFKE EENQKCKKSG GKYSLNSGAV KNPKTQLGQK PFTCSVCGKG FSQSANLVVH QRIHTGEKPF ECHECGKAFI QSANLVVHQR IHTGQKPYVC SKCGKAFTQS SNLTVHQKIH SLEKTFKCNE CEKAFSYSSQ LARHQKVHIT EKCYESNECG KTFTRSSNLI VHQRIHTGEK PFACNDCGKA FTQSANLIVH QRSHTGEKPY ECKEKGKAFS CFSHLIVHQR IHTAEKPYDC SEC GKAFSQL SCLIVHQRIH SGDLPYVCNE CGKAFTCSSY LLIHQRIHNG EKPYTCNECG KAFRQRSSLT VHQRTHTGEK PYECEKCGAA FISNSHLMRH HRTHLVE
Source	Mammalian cell
Target Names	ZNF35
Protein Names	Recommended name: Zinc finger protein 35 Alternative name(s): Zinc finger protein HF.10
Expression Region	1-527
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.