



# Recombinant Escherichia coli Ribosomal protein L11 methyltransferase (prmA)

<b>Product Code</b>	CSB-MP457802ENT
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B1XHM8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12 / DH10B)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPWIQLKLNT TGANAEDLSD ALMEAGAVSI TFQDTHDTPV FEPLPGETRL WGDTDVIGLF DAETDMNDVV AILENHPLLG AGFAHKIEQL EDKDWEREWM DNFHPMRFGE RLWICPSWRD VPDENAVNVM LDPGLAFGTG THPTTSLCLQ WLDSLDTGK TVIDFGCGSG ILAIAALKLG AAKAIGIDID PQAIQASRDN AERNGVSDRL ELYLPKDQPE EMKADV VVAN ILAGPLRELA PLISVLPVSG GLLGLSGILA SQAESVCEAY ADSFALDPVV EKEEWCRITG RKN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	prmA
<b>Protein Names</b>	Recommended name: Ribosomal protein L11 methyltransferase Short name= L11 Mtase EC= 2.1.1.-
<b>Expression Region</b>	1-293
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.