



Recombinant Thermoproteus neutrophilus Signal recognition 54 kDa protein (srp54)

Product Code	CSB-EP457878TNE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1Y9L4
Product Type	Recombinant Protein
Immunogen Species	Pyrobaculum neutrophilum (strain DSM 2338 / JCM 9278 / V24Sta) (Thermoproteus neutrophilus)
Purity	≥85% (SDS-PAGE)
Sequence	MKALSEVFGK LVEKIRGVSY IDEATLQDLS REIQRALLKA DVPLDMVKAF TDAAVKRIKE EKPPAGIPPR EYLLYVLYEE LVKLLGGEQP PEFKPTKKPY VLLLLGVEGS GKTTTAAKLA KLLVKRGYKV GLVETDTPVRP AAFDQLRQLA EKIGVPFYGE RDGKDAVEIA RRGVQNLKLN DVVIVDTAGR HRNEEALLQE VKAIYEAVSP DEVVLVIDAT VGKLAQAQAE AFMRYLPIHS VIITKMDSTA RGGGALAAVA KTGAKVKFIG VGEDVDELEM FIPRKFFVARV LGMGDLDALV ERIKAVFEEE QVIQEIESGK LDLLTFKKQI DGLLKLGPLS KVFQMLPGGL AAKISEEQIE LSQKNLKKWR AILSSMTVEE LKNPELLNAS RIRRIALGAG VTPRDVKEML TVYENLKRMS KTLKRQLRMK MPR
Source	E.coli
Target Names	srp54
Protein Names	Recommended name: Signal recognition 54 kDa protein Alternative name(s): SRP54
Expression Region	1-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.