



Recombinant Thermoproteus neutrophilus 50S ribosomal protein L22P (rpl22p)

Product Code	CSB-YP457890TNE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B1YD24
Product Type	Recombinant Protein
Immunogen Species	Pyrobaculum neutrophilum (strain DSM 2338 / JCM 9278 / V24Sta) (Thermoproteus neutrophilus)
Purity	>85% (SDS-PAGE)
Sequence	MPRHQNYSL S DEAAVQLVFR KYGVKVS AEQ IAKAYAPEQK MSWKKSVEVA RFIEG MTLRQ AKSWLEDVVK MKRPIPIRTF KKKQAHHATP WSGWPVAKWP VKVARRFLDL LENLENNARF RGLDVERVVI VHAAAHKGYR IHNIMPRAFG RATRFDEQTV NVEVIAAELP QQVVPKRYRL NLV KR
Source	Yeast
Target Names	rpl22
Protein Names	Recommended name: 50S ribosomal protein L22P
Expression Region	1-185
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.