



Recombinant *Drosophila erecta* Ribosome biogenesis protein WDR12 homolog (GG10358)

Product Code	CSB-MP459253DLK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3N534
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDVENGEGQV QVHLKTKQEH YAVPDVPYAI DGTVTTVELN TFDVALLRQK DGSSDTEFDF LVFDEYLRGR LCDHLREKAI SFEDAIEIY VERFPAPEPQ DCLLHDDWVS AVKASGKWIL SGCYDNTLNL WTNKKGHILT ISGHTAPIKA VDWISLDEET GRFVSTSQDQ TAMLWQWNVG SNSVECVSVC KGERGVDSV SVSPDGLRFA TGSWDTMLKV WSAEQEDAAE GSSKRMKESG VRTPKITLQG HRESVSAVQW MDASTLLTGS WDHTLKVWDL SLEGIKTEIS TNKSIFDASY SKLNRLILTA SADKNLRLYD PRTNQGSVVR NTYLGHNAWV QTVMWSTTEE FLFVSGAYDN QNKLWDCRSP KAPLYDLLGH GEKVLIDIDWS NPKYIVSGGV DNTVRVFKSR KALAEDEAK
Source	Mammalian cell
Target Names	GG10358
Protein Names	Recommended name: Ribosome biogenesis protein WDR12 homolog
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.