



Recombinant *Drosophila sechellia* Probable cytosolic Fe-S cluster assembly factor GM20417 (GM20417)

| | |
|--------------------------|--|
| Product Code | CSB-BP459611DMG |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B4IMH3 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Drosophila sechellia</i> (Fruit fly) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSRWSTALQL TDIDDFITPS QICIKPVQID KARSKTGAKI KIKGDSCFEE SESGNLKLNK VDISLQDCLA CSGCITSAEE VLITQQSQEE LLKILQENSK NKASEDWDNV RTIVITLATQ PLLSLAHRYQ IGVEDAARHL NGYFRSLGAD YVLSTKVADD IALLECRQEF VDRYRENEENL TMLSSSCPGW VCYAEKTHGN FILPYVSTTR SPQQIMGVLV KQILADKINV PASRIYHVTM MPCYDKKLEA SREDFFSKAN NSRDVDCVIT SVEVEQLLSE AQRPLSQYDL FDLWDPWSNV RPEFMVWAHE KTQSGGYAEH IFKFAAKHIF NEDLTTELEF KQLKNRDFRE IILKQNGKTV LKFAIANGFR NIQNLVQKLG REKLSNYHFV EVMACPSGCI NGGAQIRPTT GQHVRELTRK LEELYHNLPR SEPENSLTKH IYNDFLDGFQ TDKSYEMLHT RYHDVVSDLS ISLNINW |
| Source | Baculovirus |
| Target Names | GM20417 |
| Protein Names | Recommended name: Probable cytosolic Fe-S cluster assembly factor GM20417 |
| Expression Region | 1-477 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.