



Recombinant Human Endoplasmin (HSP90B1)

Product Code	CSB-BP010814HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P14625
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	DDEVVDV DGT VEEDLGKSRE GSRTDDEVVQ REEEAIQLDG LNASQIRELR EKSEKFAFQA EVNRMMKLII NSLYKNKEIF LRELISNASD ALDKIRLISL TDENALSGNE ELTVKIKCDK EKNLLHVTD T GVGMTREELV KNLGTIAKSG TSEFLNKMTE AQEDGQSTSE LIGQFGVGFY SAFLVADKVI VTSKHNNDTQ HIWESDSNEF SVIADPRGNT LGRGTTITLV LKEEASDYLE LDTIKNLVKK YSQFINFPIY VWSSKTETVE EPMEEEEAAK EEKEESDDEA AVEEEEEEEK PKTKKVEKTV WDWELMNDIK PIWQRPSKEV EEDEYKAFYK SFSKESDDPM AYIHFTAEGE VTFKSILFVP TSAPRGLFDE YGSKKSDYIK LYVRRVFITD DFHDMMPKYL NFVKGVVDS D DLPLNVSRET LQQHKLLKVI RKKLVRKTL D MIKKIADDKY NDTFWKEFGT NIKLGVIEDH SNRTRLAKLL RFQSSHHPTD ITSLDQYVER MKEKQDKIYF MAGSSRKEAE SSPFVERLLK KGYEVIYLTE PVDEYCIQAL PFDGKRFQN VAKEGVKFDE SEKTKESREA VEKEFEPLL N WMKDKALKDK IEKAVVSQRL TESPALVAS QYGWSGNMER IMKAQAYQTG KDISTNYYAS QKKTFEINPR HPLIRDMLRR IKEDEDDKTV LDLAVVLFET ATLRSGYLLP DTKAYGDRIE RMLRSLNID PDAKVEEPEE EPEETAEDT TEDTEQDEDE EMDVGTDEEE ETAKESTAEK DEL
Source	Baculovirus
Target Names	HSP90B1
Protein Names	Recommended name: Endoplasmin Alternative name(s): 94 kDa glucose-regulated protein Short name= GRP-94 Heat shock protein 90 kDa beta member 1 Tumor rejection antigen 1 gp96 homolog
Expression Region	22-803
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.