



# Recombinant Human Endoplasmin (HSP90B1)

<b>Product Code</b>	CSB-MP010814HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P14625
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	DDEVVDVGT VEEDLGKSRE GSRTDDEVVQ REEEAIQLDG LNASQIRELR EKSEKFAFQA EVNRMMKLII NSLYKNKEIF LRELISNASD ALDKIRLISL TDENALSGNE ELTVKIKCDK EKNLLHVTDG GVGMTREELV KNLGTIAKSG TSEFLNKMTE AQEDGQSTSE LIGQFGVGFY SAFLVADKVI VTSKHNNDTQ HIWESDSNEF SVIADPRGNT LGRGTTITLV LKEEASDYLE LDTIKNLVKK YSQFINFPIY VWSSKTETVE EPMEEEEAAK EEKEESDDEA AVEEEEEEEK PKTKKVEKTV WDWELMNDIK PIWQRPSKEV EEDEYKAFYK SFSKESDDPM AYIHFTAEGE VTFKSILFVP TSAPRGLFDE YGSKKSDYIK LYVRRVFITD DFHDMMPKYL NFVKGVVDS DPLPLNVSRET LQQHKLLKVI RKKLVRKTL D MIKKIADDKY NDTFWKEFGT NIKLGVIEDH SNRTRLAKLL RFQSSHHPTD ITSLDQYVER MKEKQDKIYF MAGSSRKEAE SSPFVERLLK KGYEVIYLTE PVDEYCIQAL PFDGKRFQN VAKEGVKFDE SEKTKESREA VEKEFEPLL N WMKDKALKDK IEKAVVSQRL TESPALVAS QYGWSGNMER IMKAQAYQTG KDISTNYYAS QKKTFEINPR HPLIRDMLRR IKEDEDDKTV LDLAVVLFET ATLRSGYLLP DTKAYGDRIE RMLRSLNID PDAKVEEPEE EPEETAEDT TEDTEQDEDE EMDVGTDEEE ETAKESTAEK DEL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	HSP90B1
<b>Protein Names</b>	Recommended name: Endoplasmin Alternative name(s): 94 kDa glucose-regulated protein Short name= GRP-94 Heat shock protein 90 kDa beta member 1 Tumor rejection antigen 1 gp96 homolog
<b>Expression Region</b>	22-803
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.