



Recombinant *Drosophila ananassae* Cytosolic Fe-S cluster assembly factor NUBP1 homolog (GF20875)

Product Code	CSB-BP461196DKY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3MU92
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MEAPPEHCPCG VESEQAGKVS ACAGCPNQGI CSDPNRKVED PGKALVAESL KDVKNKLLIL SGKGGVGKST VTSLTRYLA RSCPNSNFGV LDIDICGPSQ PRLGALGEN VHQSGSGWSP VGIDDNVCLM SIGFLLGSVD DAIWRGPKK NGMIRQFLSE VDWGNLDDLL LDTPPGTSDE HLSVVSYLKD DTQPESLRAV IVTTPQEVAL LDVRKEINFC KKQNIPIVGV IENMSSFRCG NCGNSSEIFP AKTGGAAMC AEMEVPLLGS LPLDPQIAKA CDSGDDITEI KNSTTEALDG ICSKIMSSFS
Source	Baculovirus
Target Names	GF20875
Protein Names	Recommended name: Cytosolic Fe-S cluster assembly factor NUBP1 homolog
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.