



Recombinant *Drosophila yakuba* Protein phosphatase PTC7 homolog fig (fig)

Product Code	CSB-BP461612DMR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4PPK3
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila yakuba</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MITHLKNWPR LLNRFVLQLK NARHSIHQFT HLAGRLQRPP KSGKSSRD PY LVTAVQGRSK KPRYPGERAN QRFGEDSWFV SSTPLAEVMG VADGVGGWRD VGVDAGRFAK ELMTCCSGQT QRSGFDGRSA RNLLIAGFQE LTHREQPVVG SSTACLATMH RRDCILYTAN LGDSGFLVVR NGRVLHRSVE QTHDFNTPYQ LTVPPADRQD CYYCDKPEMA VSTRHSLLP G DLVLLATDGL FDNMPESMLL KILNGLKERG ERDLLQGASQ VVEKARELSL NATFQSPFAI KARQHNPYS GGGKPDDITL ILASVEVPRA
Source	Baculovirus
Target Names	fig
Protein Names	Recommended name: Protein phosphatase PTC7 homolog fig Alternative name(s): Fos intronic gene protein EC= 3.1.3.16
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.