



Recombinant *Saccharomyces cerevisiae* 6-phosphogluconolactonase 3 (SOL3)

Product Code	CSB-EP463129SVP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3LSS7
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain RM11-1a) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MVTGVFSEER ASLTHQLGEF IVKKQDEALQ KKSDFKVSVS GGSLIDALYE SLVADESLSS RVQWSKWQIY FSDERIVPLT DADSNYGAFK RAVLDKLPST SQPNVYPMDE SLIGSDAESN NKIAAEYERI VPQVLDLVLL GCGPDGHTCS LFPGETHRYL LNETTKRVAW CHDSPKPPSD RITFTLPVLK DAKALCFVAE GSSKQNMHE IFDLKNDQLP TALVNKLFGE KTSWVFVNEEA FGKVQTKTF
Source	<i>E.coli</i>
Target Names	SOL3
Protein Names	Recommended name: 6-phosphogluconolactonase 3 Short name= 6PGL EC= 3.1.1.31
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.