



Recombinant *Drosophila ananassae* Eukaryotic translation initiation factor 3 subunit F-1 (eIF3-S5-1)

| | |
|--------------------------|---|
| Product Code | CSB-YP463136DKY |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B3M123 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Drosophila ananassae</i> (Fruit fly) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSALNLTVRV HPVVLFQVVD AYERRNADSH RVIGTLLGSV DKGVVEVTNC FCVPHKEHDD QVEAELSYAL DMYDLNRKVN PNESVVGWWA TGNDVTNHSS VIHEYAREC NNPVHLTVDT SLQGGRMGLR AYVCIQLGVP GGKSGCMFTP IPVELTSYEP ETFGLKLLQK TVGVSPAHRP KTVPPMLDLA QISEASQKLQ QLLDLILKYV DDVIAHKVTP DNAVGRQLLD LIHSVPHMTH EQFTQMFNAN VRNLLLVTIL SQLIKTQLQL NEKLTFLPSA |
| Source | Yeast |
| Target Names | eIF3f1 |
| Protein Names | Recommended name: Eukaryotic translation initiation factor 3 subunit F-1 Short name= eIF3f-1 Alternative name(s): Eukaryotic translation initiation factor 3 subunit 5-1 |
| Expression Region | 1-280 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |