



Recombinant *Drosophila ananassae* Ribosome biogenesis protein WDR12 homolog (GF14067)

Product Code	CSB-EP463140DKY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3MJV8
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MDVENGEGQV QVHLKTKQEH YAVPDVPYAI DGTVTTAELN Tfvntllqqk DGSPVDFDFL VFDEYLRGRL CDHLREKAIS FEDAIEIEYV ERFPAPEPQD CLLHDDWVSA VKANGKWILT GCYDNTLNIW TNKGKHILTI SGHTAPIKAV DWISLDEESG RFVSSSQDQT AMLWQWNVGS NTAECVSVCK GHERGVDSVS VSPDGQRFAT GSWDTMLKVW SAELDDAGEG TSKRMKESGV RTPKMTLQGH RESISAVQWM DASTLLTGSW DHTLKVWDLS LEGIKTEIST NKSIFDASYS NLNRLIVTAS ADKNLRLYDA RTNQGSSVVRN TYLGHNAWVQ AVMWSTTEEF Lfvsgaydnq NKLWDCRSPK APLYDLLGHG EKVLdidwtn PKYIVSGGAD NSRVVFKSRK AVVESMDTK
Source	E.coli
Target Names	GF14067
Protein Names	Recommended name: Ribosome biogenesis protein WDR12 homolog
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.