



Recombinant *Drosophila ananassae* Anamorsin homolog (I (2)35Bg)

Product Code	CSB-MP463141DKY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B3MLC5
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila ananassae</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MEQFKGLQKT LYIWTDSDL DKRVEKVKSA TGGEVALENV HRLSFSSYAN SSFDLIVIEC AQLTDNYVKL LHMLKPSGKL HLVAFIGPAG SLLQEVKLSG FINCREDADD ALTAEKPGYE TGSSARLSFA KKDASALNVW KISGDDEELI DEEDLLDEED KQKPDPAGLR VCSTTGKRKA CKNCSCGLAE ELESEKQSKA ATENAKSSCG NCYLGDAFRC STCPYLGMPA FKPGEKVQLA DNLLKSDI
Source	Mammalian cell
Target Names	CIAPIN1
Protein Names	Recommended name: Anamorsin homolog Alternative name(s): Fe-S cluster assembly protein DRE2 homolog
Expression Region	1-248
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.