



Recombinant Drosophila sechellia SAGA-associated factor 11 homolog (Sgf11)

Product Code	CSB-BP463508DMG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B4IFU5
Product Type	Recombinant Protein
Immunogen Species	Drosophila sechellia (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSAANMPTTT GAQGSGNQVP RTSTTIVNHR ELIKEPKNLD EAANYLFQTL LDDAVVGIFN ETHHLRKSGN LAALDGVPED STYRMCEMPN LDIFGISTAK KPMDCPCPNC DRLVAAARFA PHLEKCMGMG RISSRIASRR LATKEGATSA HLHSSGNTGG TDEDDVDWS SDKRRKKSQ NSRNGSKKN NGKTF
Source	Baculovirus
Target Names	Sgf11
Protein Names	Recommended name: SAGA-associated factor 11 homolog
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.